

Companion® 31 & 41

Liquid Oxygen Reservoirs

The One that Set the Standards

Companion® Liquid Oxygen Reservoirs offer high quality and reliability. Companion 31 and 41 models have capacities of 31 and 41 liquid liters of oxygen, respectively. The proven vessel design of these reservoirs is rugged and easy to use. Low maintenance and long life make them cost-effective to own. The Companion 31 and 41 reservoirs are compatible with all Companion Portable Liquid Oxygen Units.

Features

- Low maintenance; no batteries required
- Enhanced insulation system decreases evaporative oxygen losses
- Built-in flow control valve provides 11 standard settings from 1/4 to 6 LPM
- Mechanical contents indicators give reliable reading of remaining oxygen levels



Call 1.800.482.2473 to order!
www.CAIREmedical.com

Companion® 31 & 41

Liquid Oxygen Reservoirs



Product Specifications

	Companion 31	Companion 41
LOX capacity	73 lb (33,11 kg)	98 lb (44,45 kg)
Gaseous equivalent capacity	24.950 L	33,500 L
Weight, empty	51 lb (23,13 kg)	60 lb (27,21 kg)
Weight, filled	124 lb (56,25 kg)	158 lb (71,67 kg)
Height	33" (838 mm)	38.5" (978 mm)
Diameter	14.2" (361 mm)	14.2" (361 mm)
Typical use time at 2 LPM 16 BPM demand mode	8 days 20 hrs	11 days 17 hrs
Operating pressure	19.5 psi (134 kPa)	19.5 psi (134 kPa)
Normal evaporation rate	1.9 lb/day (0,86 kg/day)	1.9 lb/day (0,86 kg/day)
Standard flow control range	Off, .25, .75, 1, 1.5, 2, 2.5, 3, 4, 5, 6	Off, .25, .75, 1, 1.5, 2, 2.5, 3, 4, 5, 6

* Standard fill times vary with the saturation and operating pressures of the source vessel. Contact Technical Service for more information.

Ordering Information*

Description	Catalog Number
Companion 31A Reservoir	B-700893-00
Companion 41A Reservoir	B-700895-00

Accessories

Roller Base	10855678
Transfer Line Assembly, 6'	B-775288-00
Transfer Line Assembly, 10'	B-775289-00
Universal Adapter Kit	B-775461-00
Pressure Gauge, tube adapter	B-775270-00

* Catalog Numbers reflect US English Part Numbers, for language-specific part numbers please contact Customer Service.

CAIRE, Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107
Ph 770-721-7759 • Toll-Free 1-800-482-2473
Fax 888-932-2473

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham RG41 2PL, United Kingdom
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245