

Liquid Oxygen Duration Chart

Portables

	Off	0.25	0.5	0.75	1	1.5	2	2.5	3	4	5	6	8	10	15
C1000	54.5	30.8	21.5	16.5	13.4	9.7	7.7	6.3	5.4	4.1	3.3	2.8			
C1000T	54.5		21.5		13.4	9.7	7.7		5.4	4.1	3.3	2.8	2.1	1.7	1.2

* Portable time is listed in hours

Bases

	Off	0.25	0.5	0.75	1	1.5	2	2.5	3	4	5	6	8	10	15
Companion 41	51.8	51.8	44.2	34.5	25.2	17.1	12.7	10	8.3	6.2	5	4.1	2.9	2.7	
MVE 45	74.6	74.6	58	37.8	28.1	18.8	14.1	11.3	9.4	7	5.6	4.7	3.5	2.8	1.9
MVE 60	90	90	74.7	48.7	36.2	24.2	18.1	14.5	12.1	9.1	7.2	6	4.5	3.6	2.4

* Base time is listed in days