





DELUXE CPAP MASK SIZE GUIDE

877-578-6738
www.SunsetHCS.com


LARGE




MEDIUM



NASAL MASK
CM006




Mouth should be slightly open,
as if breathing.



FULL FACE MASK
CM007

IMPORTANT Print without page scaling (print at 100% scale). Turn off any auto scaling or shrink to fit options.
Printed area should measure 8 inches by 10.5 inches. Cut inside blue and red guidelines. Place sheet over face.

LARGE



MEDIUM



SMALL

